
BioRipar[®]


CRPATCH[®] - Cystocele Repair Patch[®]

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Pelvic organic support, described by J. DeLancey, (J.O.L. DELANCEY 1994) identifies in the level II the anchorage of the medial portion of the vaginal walls through the **pubo-cervical fascia**, which is fixed to the **tendinous arch of the endo-pelvic fascia**. Therefore, the anatomic integrity is the main mandatory condition to the functional restore of pelvic organs.

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
Standard pelvic floor surgery, main phases:

- I. Anterior Vagina wall Incision «reverse T» form the vesical neck to the portio, bilaterally, or to the vaginal cupola, in case of previous hysterectomy
- II. Pubic-vesical fascia separation or "Alban's fascia" from the vagina wall to the ischium-pubic branch
- III. Pubic-vesical fascia duplication thorough monofilament long term absorbable suture Cal. 2/0 (Assufil Monofilament® FU395M) separate stich technique.
- IV. Vaginal wall recontation
- V. Vaginal Breach closure thought continuous stitch technique thorough barbed long term absorbable suture (Filbloc® FQ07GHHAD).

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CRPatch® Cystocele Repair Patch® is made of pure collagen of bovine pericardium tissue, **non cross-linked**, which favorites tissue remodeling process through the physical organic interaction with the vaginal wall. With its anti-infectious properties, it guarantees an absence of complications generally found in the use of synthetic prostheses (*safety*)*.

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
The integration in the host tissue three months after implantation, shows an advanced integration, with evident deposit of neo-collagen and angiogenesis. The process will finish with a physiological regeneration of a new functional structure with increased tensile strength, useful to repair the anatomic deficit of cystocele, the compartment with the highest relapse rate in pelvic surgery.

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CRPatch® Cystocele Repair Patch® has a specific design developed to optimize the surgical strategies oriented to the functional anatomy in terms of resolution of the specific defect. The shape has been designed to recreate the anatomy of DeLancey's level II, reinforcing specifically the damaged tissues, fixing the patch's wings to the tendinous arch of the endo-pelvic fascia. The sub-urethral positioning and the **CRPatch®**'s medial body perform the repair of the central defect of the cystocele, the lateral wings are used to solve the lateral defect by functional anchorage.


The specific procedure is the same compared to the standard surgery for the first 3 steps, then, the paravesical spaces are prepared, for a depth of about 3 cm, towards the tendinous arch, allocating the lateral wings in this site, immediately after the sub-urethral positioning. The fixation to the pubo-cervical fascia is performed with Cal. 2/0 monofilament long term absorbable suture (Assufil Monofilament® FU395M) with separate stitch technique.


Bibliografia - Bibliography

- Olsen AL, Smith VJ, Bergstrom JO, Colling JC, Clark AL. Obstet Gynecol. 1997 Apr;89(4):501-6
- Birch C, Fynes MM. Curr Opin Obstet Gynecol. 2002 Oct; 14(5): 527-35.
- Bump RC, Mattiasson A, Bø K, Brubaker LP, DeLancey JO, Klarskov P, Shull BL, Smith AR. Am J Obstet Gynecol. 1996 Jul;175(1):10-7.
- Lien YS1, Chen GD, Ng SC. Int J Gynaecol Obstet. 2012 Nov;119(2):185-8
- Maher C, Feiner B, Baessler K, Adams EJ, Hagen S, Glazener CM. Cochrane Database Syst Rev. 2010 Apr 14;(4)
- FDA CDRH Urogynecologic Surgical Mesh: Update on the Safety and Effectiveness of Transvaginal Placement for Pelvic Organ Prolapse July 2011.7.
- Konstantinovic ML, Lagae P, Zheng F, Verbeken EK, De Ridder D, Deprest JA. BJOG. 2005 Nov;112(11):1554-60.
- Rufus Cartwright, MA Systematic review and metaanalysis of genetic association studies of urinary symptoms and prolapse in women American Journal of Obstetrics & Gynecology FEBRUARY 2015
- NHS England Mesh Oversight Group Report 25 July 2017
- Carl J Heneghan Trials of transvaginal mesh devices for pelvic organ prolapse: a systematic database review of the US FDA approval process BMJ Open 2017

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© 2019 Assut Europe S.p.A.all rights reserved

CE 0373

UNI EN ISO 9001:2015 - UNI EN ISO 13485:2016 - FDA 9617547

This document is intended solely for the use of healthcare professionals.

Please refer always to the instruction for use.

Products may not be available in all markets.

Please contact your Assut Europe representative if you have questions about the availability of products in your area